

# Product Information

## **MemDX™ Membrane Protein Human BCAP29 (B cell receptor associated protein 29)**

### **Expressed *in vitro* E.coli expression system, Full Length**

Cat. No.: **MPX1459K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human BCAP29 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### **Product Specifications**

#### **Host Species**

Human

#### **Target Protein**

BCAP29

#### **Protein Length**

Full Length

#### **Protein Class**

Transport

#### **TMD**

3

#### **Sequence**

MTLQWAAVATFLYAEIGLILIFCLPFIPPQRWQKIFSFNVWGKIATFWNKAFLTIILLIVLFLDAVREVRKYSSVHTIEKSSTSRPDAYEH

### **Product Description**

#### **Expression Systems**

*in vitro* E.coli expression system

#### **Tag**

10xHis tag at the N-terminus

#### **Protein Format**

Soluble

#### **Form**

Liquid or Lyophilized powder

#### **Buffer**

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

**Storage**

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

**Target****Target Protein**

BCAP29

**Full Name**

B cell receptor associated protein 29

**Alternative Names**

BCAP29; BAP29; B-cell receptor-associated protein 29; BCR-associated protein 29; B cell receptor associated protein 29

**Gene ID**

[55973](#)

**UniProt ID**

[Q9UHQ4](#)