

Product Information

MemDX™ Membrane Protein Human DNAJC4 (DnaJ heat shock protein family (Hsp40) member C4) Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX1461K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human DNAJC4 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

DNAJC4

Protein Length

Full Length

Protein Class

Receptor

TMD

1

Sequence

MPPLLPLRLCRLWPRNPPSRLLGAAAGQRSRPSTYYELLGVHPGASTE EVKRAFFSKSKELHPDRDPGNPSLHSRFVELSEAYRVI

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

DNAJC4

Full Name

DnaJ heat shock protein family (Hsp40) member C4

Alternative Names

DNAJC4; HSPF2; MCG18; DANJC4; dnaJ homolog subfamily C member 4; DnaJ (Hsp40) homolog, subfamily C, member 4; dnaJ-like protein HSPF2; heat shock 40kD protein 2; multiple endocrine neoplasia type 1 candidate protein number 18; DnaJ heat shock protein family (Hsp40) member C4

Gene ID

[3338](#)

UniProt ID

[Q9NNZ3](#)