

Product Information

MemDX™ Membrane Protein Human TMEM174 (Transmembrane protein 174) Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX1465K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human TMEM174 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TMEM174

Protein Length

Full Length

Protein Class

Receptor

TMD

2

Sequence

MEQGSGRLEDFPVNVFSVTPYTPSTADIQVSDDDKAGATLLFSGIFLGLVGITFTVMGWIKYQGVSHFEWTQLLGPVLLSVGVTFIL

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

TMEM174

Full Name

Transmembrane protein 174

Alternative Names

TMEM174; FLJ31268, MGC13034; Transmembrane protein 174

Gene ID

[134288](#)

UniProt ID

[Q8WUU8](#)