

Product Information

MemDX™ Membrane Protein Mouse Rnf182 (Ring finger protein 182) Expressed *in vitro* *E.coli* expression system, Full Length

Cat. No.: **MPX1474K**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse Rnf182 membrane protein expressed *in vitro E.coli* expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Mouse

Target Protein

Rnf182

Protein Length

Full Length

Protein Class

Transferase

TMD

2

Sequence

MASQPLEEPAESQASDELECKICYNRYNLKQRKPKVLECCHRVCAKCLYKIIDFGDSPQGVIVCPFCRFETCLPDDEVSSLPDDNNI

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

Rnf182

Full Name

Ring finger protein 182

Alternative Names

Rnf182; C630023L15Rik; E3 ubiquitin-protein ligase RNF182; RING-type E3 ubiquitin transferase RNF182; Ring finger protein 182

Gene ID

[328234](#)

UniProt ID

[Q8C432](#)