

Product Information

MemDX™ Membrane Protein Mouse Mettl23 (Methyltransferase like 23) Expressed *in vitro* *E.coli* expression system, Full Length

Cat. No.: **MPX1495K**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse Mettl23 membrane protein expressed *in vitro E.coli* expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Mouse

Target Protein

Mettl23

Protein Length

Full Length

Protein Class

Transferase

TMD

1

Sequence

MDSVRPRAPWAPPPDPASLDSPTCEPGLMAGTQLFRFREEPVPGGNRAVLEVRVPQVLHVQYGMVWPCAVVLAQYLWFHRRS

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

Mettl23

Full Name

Methyltransferase like 23

Alternative Names

Mettl23; 1110005A03Rik; 1500035B17Rik; 4933424L15Rik; methyltransferase-like protein 23; Methyltransferase like 23

Gene ID

[74319](#)

UniProt ID

[A2AA28](#)