

Product Information

MemDX™ Membrane Protein Human REEP4 (Receptor accessory protein 4) Expressed *in vitro* *E.coli* expression system, Full Length

Cat. No.: **MPX1510K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human REEP4 membrane protein expressed *in vitro* *E.coli* expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

REEP4

Protein Length

Full Length

Protein Class

Receptor

TMD

2

Sequence

MVSWMICRLVVLVFGMLCPAYASYKAVKTKNIREYVRWMMYWIVFALFMAAEIVTDIFISWFPFYIEIKMAFVLWLLSPYTKGASLLY

Product Description

Expression Systems

in vitro *E.coli* expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

REEP4

Full Name

Receptor accessory protein 4

Alternative Names

REEP4; PP432; Yip2c; C8orf20; receptor expression-enhancing protein 4; Receptor accessory protein 4

Gene ID

[80346](#)

UniProt ID

[Q9H6H4](#)