

Product Information

MemDX™ Membrane Protein Mouse Zdhhc22 (Zinc finger, DHHC-type containing 22)

Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX1525K**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse Zdhhc22 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Mouse

Target Protein

Zdhhc22

Protein Length

Full Length

Protein Class

Transferase

TMD

5

Sequence

MLALRLLNVVAPAYFLCISLVTFVLQLFLFLPSMREDPTATPLFSPAHLHGALFLFLSANALGNYVLVIQNSPDDLGTCCQGTMSQRPG

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

Zdhhc22

Full Name

Zinc finger, DHHC-type containing 22

Alternative Names

Zdhhc22; Gm262; palmitoyltransferase ZDHHC22; DHHC-22; putative palmitoyltransferase ZDHHC22; zinc finger DHHC domain-containing protein 22; Zinc finger, DHHC-type containing 22

Gene ID

[238331](#)

UniProt ID

[A0PK84](#)