

Product Information

MemDX™ Membrane Protein Human ITM2C (Integral membrane protein 2C) Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX1533K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human ITM2C membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

ITM2C

Protein Length

Full Length

Protein Class

Receptor

TMD

1

Sequence

MVKISFQPAVAGIKGDKADKASAPAPASATEILLTPAREEQPPQHRSKRGGSVGGVCYLMSGMVLLMGLVFASVYIYRYFFLAC

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

ITM2C

Full Name

Integral membrane protein 2C

Alternative Names

ITM2C; E25; BRI3; E25C; ITM3; BRICD2C; BRICHOS domain containing 2C; cerebral protein 14; integral membrane protein 3; transmembrane protein BRI3; Integral membrane protein 2C

Gene ID

[81618](#)

UniProt ID

[Q9NQX7](#)