

Product Information

MemDX™ Membrane Protein Human RNF148 (Ring finger protein 148) Expressed *in vitro* *E.coli* expression system, Full Length of Mature Protein

Cat. No.: **MPX1539K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human RNF148 membrane protein expressed *in vitro E.coli* expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

RNF148

Protein Length

Full Length of Mature Protein

Protein Class

Receptor

TMD

1

Sequence

KAIWTAHLNITFQVGNEITSELGESGVFGNHSPLERVSGVVALPEGWNQNACHPLTNFSRPKQADSWLALIERGGCTFTHKINVAAR

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

RNF148

Full Name

Ring finger protein 148

Alternative Names

RNF148; RING finger protein 148; Ring finger protein 148

Gene ID

[378925](#)

UniProt ID

[Q8N7C7](#)