

Product Information

MemDX™ Membrane Protein Human ABHD14A (Abhydrolase domain containing 14A)

Expressed *in vitro* *E.coli* expression system, Full Length

Cat. No.: **MPX1542K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human ABHD14A membrane protein expressed *in vitro* *E.coli* expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

ABHD14A

Protein Length

Full Length

Protein Class

Receptor

TMD

1

Sequence

MVGALCGCWFRLLGGARPLIPLGPTVVQTSMSRSQVALLGLSLLLMLLLYVGLPGPPEQTSCLWGDPNVTVLAGLTPGNSPIFYREV

Product Description

Expression Systems

in vitro *E.coli* expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

ABHD14A

Full Name

Abhydrolase domain containing 14A

Alternative Names

ABHD14A; DORZ1; protein ABHD14A; abhydrolase domain-containing protein 14A; alpha/beta hydrolase domain-containing protein 14A; Abhydrolase domain containing 14A

Gene ID

[25864](#)

UniProt ID

[Q9BUJ0](#)