

Product Information

MemDX™ Membrane Protein Human FAM210A (Family with sequence similarity 210 member

A) Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX1546K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human FAM210A membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

FAM210A

Protein Length

Full Length

Protein Class

Receptor

TMD

1

Sequence

MQWNVPRTVSRLARRTCLEPHNAGLFGHCQNVKGPLLLYNAESKVVLVQGPQKQWLHLSAAQCVAKERRPLDAHPPQPGVLRH

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

FAM210A

Full Name

Family with sequence similarity 210 member A

Alternative Names

FAM210A; HsT2329; C18orf19; protein FAM210A; uncharacterized protein C18orf19; Family with sequence similarity 210 member A

Gene ID

[125228](#)

UniProt ID

[Q96ND0](#)