

Product Information

MemDX™ Membrane Protein Mouse Ergic3 (ERGIC and golgi 3) Expressed *in vitro* *E.coli* expression system, Full Length

Cat. No.: **MPX1758K**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse Ergic3 membrane protein expressed *in vitro* *E.coli* expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Mouse

Target Protein

Ergic3

Protein Length

Full Length

Protein Class

Transport

TMD

2

Sequence

MEALGKLKQFDAYPKTLED FRVKTCGGATVTIVSGLLM LLLFLSELQYYLTTEVHPELYVDKSRGDKLKINIDVLFPHMPCAYLSIDAM

Product Description

Expression Systems

in vitro *E.coli* expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

Ergic3

Full Name

ERGIC and golgi 3

Alternative Names

Ergic3; CGI-5; D2Ucl; CGI-54; NY-BR-; Sdbcag; D2Ucla1; AV318804; NY-BR-84; Sdbcag84; 2310015B14Rik; endoplasmic reticulum-Golgi intermediate compartment protein 3; serologically defined breast cancer antigen 84; serologically defined breast cancer antigen NY-BR-84 homolog; ERGIC and golgi 3

Gene ID

[66366](#)

UniProt ID

[Q9CQE7](#)