

Product Information

MemDX™ Membrane Protein Mouse Kcnj3 (Potassium inwardly-rectifying channel, subfamily J, member 3) Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX1865K**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse Kcnj3 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Mouse

Target Protein

Kcnj3

Protein Length

Full Length

Protein Class

Ion channel, Transport

TMD

2

Sequence

MSALRRKFGDDYQVVTSSSSGSLQPQGPQGQPQQQLVPKKKRQRFVDKNGRCNVQHGNLGSETSRYLSDLFTTLVDLKWRWN

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

Kcnj3

Full Name

Potassium inwardly-rectifying channel, subfamily J, member 3

Alternative Names

Kcnj3; GIR; Kcnf; Kir3; GIRK1; Kcnf3; Kcnj1; GIRK-1; Kir3.1; G protein-activated inward rectifier potassium channel 1; G protein-gated K⁺ channel subunit 1 GIRK1; inward rectifier K(+) channel Kir3.1; potassium channel, inwardly rectifying subfamily J member 3; Potassium inwardly-rectifying channel, subfamily J, member 3

Gene ID

[16519](#)

UniProt ID

[P63250](#)