

Product Information

MemDX™ Membrane Protein Human RETREG2 (Reticulophagy regulator family member 2)

Expressed *in vitro* *E.coli* expression system, Full Length

Cat. No.: **MPX1904K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human RETREG2 membrane protein expressed *in vitro* *E.coli* expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

RETREG2

Protein Length

Full Length

Protein Class

Receptor

TMD

3

Sequence

MASGGGGGNTGAGGGPGMGLSLGLGLSLGMSEATSEAAAAATAEAVGRLATTLWLRLRGWEAVLAAAQRLLVWEKPLHSLV

Product Description

Expression Systems

in vitro *E.coli* expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

RETREG2

Full Name

Reticulophagy regulator family member 2

Alternative Names

RETREG2; MAG-2; C2orf17; FAM134A; family with sequence similarity 134 member A; metastasis associated gene-2; protein FAM134A; Reticulophagy regulator family member 2

Gene ID

[79137](#)

UniProt ID

[Q8NC44](#)