

Product Information

MemDX™ Membrane Protein Human NEPRO (Nucleolus and neural progenitor protein)

Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX1925K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human NEPRO membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

NEPRO

Protein Length

Full Length

Protein Class

Receptor

Sequence

MMAAVPPGLEPWNRVRIPKAGNRSVTVQNPGAALDLCIAAVIKECHLVILSLKSQTLDAETDVLCAVLYSNHNRMGRHKPHLALKQ

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

NEPRO

Full Name

Nucleolus and neural progenitor protein

Alternative Names

NEPRO; ANXD3; NET17; C3orf17; protein nepro homolog; Nucleolus and neural progenitor protein

Gene ID

[25871](#)

UniProt ID

[Q6NW34](#)