

Product Information

MemDX™ Membrane Protein Human LRRC70 (Leucine rich repeat containing 70) Expressed *in vitro* E.coli expression system, Full Length of Mature Protein

Cat. No.: **MPX1944K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human LRRC70 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

LRRC70

Protein Length

Full Length of Mature Protein

Protein Class

Receptor

TMD

1

Sequence

CSSVCQLCTGRQINCRNLGLSSIPKNFPESTVFLYLTGNNISYINESELTGLHSLVALYLDNSNILYVYPKAFVQLRHLYFLFLNNNFIK

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

LRRC70

Full Name

Leucine rich repeat containing 70

Alternative Names

LRRC70; SLRN; leucine-rich repeat-containing protein 70; synleurin; Leucine rich repeat containing 70

Gene ID

[100130733](#)

UniProt ID

[Q7Z2Q7](#)