

Product Information

MemDX™ Membrane Protein Mouse Lrrc4c (Leucine rich repeat containing 4C) Expressed *in vitro* E.coli expression system, Full Length of Mature Protein

Cat. No.: **MPX1949K**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse Lrrc4c membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Mouse

Target Protein

Lrrc4c

Protein Length

Full Length of Mature Protein

Protein Class

Receptor

TMD

1

Sequence

QTCPSVCSCSNQFSKVICVRKNLREVPDGISTNTRLLNLHENQIQIIKVNFSFKHLRHLEILQLSRNHIRTIEIGAFNGLANLNTLELFDNF

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

Lrrc4c

Full Name

Leucine rich repeat containing 4C

Alternative Names

Lrrc4c; NGL-1; 6430556C10Rik; leucine-rich repeat-containing protein 4C; netrin-G1 ligand; Leucine rich repeat containing 4C

Gene ID

[241568](#)

UniProt ID

[Q8C031](#)