

Product Information

MemDX™ Membrane Protein SARS2 S Mutation (Δ H69/V70, Δ Y144, E484K, N501Y, A570D, P681H, T716I, S982A, D1118H, D614G) Protein (Surface glycoprotein) Reporter Virus

Particles

Cat. No.: **MPX2001K**

This product is for research use only and is not intended for diagnostic use.

This product is a SARS2 S Mutation (Δ H69/V70, Δ Y144, E484K, N501Y, A570D, P681H, T716I, S982A, D1118H, D614G) protein Reporter Virus Particles. Reporter Virus Particles (RVPs) are replication-incompetent pseudotyped virus particles. RVPs display antigenically correct virus protein on a heterologous virus core and carry a modified genome that expresses a convenient optical reporter gene (GFP or luciferase). The product is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

SARS2

Target Protein

S Mutation (Δ H69/V70, Δ Y144, E484K, N501Y, A570D, P681H, T716I, S982A, D1118H, D614G)

Protein Class

Host-virus interaction

TMD

1

Sequence

MFVFLVLLPLVSSQCVNLTTRTQLPPAYTNSFTRGVYYPDKVFRSSVLHS
TQDLFLPFFSNVTWFHAISGTNGTKRFDNPVLPFNDGVYFASTEKSN
IRGWIFGTTLDSTQSLIVNATNVVIVKCEFCNDPFLGVYHKNNK
SWMESEFRVYSSANNCTFEYVSQPFLMDLEGKQGNFKNLREFVFKNIDGY
FKIYSKHTPINLVRDLPPQGFSALEPLVDLPIGINITRFQTLALHRSYLT
PGDSSSGWTAGAAAYYVGYLQPRTFLLKYNENGTITDAVDCALDPLSETK
CTLKSFTVEKGIYQTSNFRVQPTESIVRFPNITNLCPFGEVFNATRFASV
YAWNKRISNCVADYSVLVNSASFSTFKCYGVSPKTLNDLCFTNVYADSF
VIRGDEVQRQIAPGQTGKIADYNYKLPDDFTGCVIAWNSNNLDSKVGGNYN
YLYRLFRKSNLKPFRDISTEIQAGSTPCNGVKGFNCYFPLQSYGFQPT
YGVGYQPYRVVLSFELLHAPATVCGPKKSTNLVKNKCVNFNFNGLTGTG
VLTESNKKFLPFQQFGRDIDDTDAVRDPQTLEILDITPCSFGGVSVITP
GTNTSNQVAVLYQGVNCTEVPVAIHADQLTPTWRVYSTGSNVFQTRAGCL
IGAELVNNSECDIPIGAGICASYQTQTNSHRRARSVASQSIAYTMSLG
AENSVAYSNNIAIPNFTISVTTEILPVSMTKTSVDCTMYICGDSTEC
NLLLQYGSFCTQLNRALTGIAVEQDKNTQEVFAQVKQIKTPPIKDFGGF
NFSQILPDPSKPSKRSFIEDLLFNKVTADAGFIKQYGDCLGDIAARDLI
CAQKFNGLTVPPLLTDEMIAQYTSALLAGTITSGWTFGAGAALQIPFAM
QMAYRFNGIGVTQNVLYENQKLIANQFNSAIGKIQDSLSTASALGKLQD
VVNQNAQALNTLVKQLSSNFGAIVSVLNDILARLDKVEAEVQIDRLITGR
LQSLQTYVTQQLIRAAEIRASANLAATKMSECVLGQSKRVDFCGKGYHLM

SFPQSAPHGVVFLHVTYVPAQEKNFTTAPAICHGDKAHFPREGVFVSNGT
HWFVTQRNFYEPQIITTHNTFVSGNCDVVIGIVNNTVYDPLQPELDSFKE
ELDKYFKNHTSPDVLGDISGINASVVNIQKEIDRLNEVAKNLNESLIDL
QELGKYEQYIKWPWYIWLGFIAGLIAIVMTIMLCCMTSCCCLKGCCSC
GSCCKFDEDDSEPVLKGVKLHYT

Product Description

Protein Format

Reporter Virus Particles

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

S Mutation (ΔH69/V70, ΔY144, E484K, N501Y, A570D, P681H, T716I, S982A, D1118H, D614G)

Full Name

Surface glycoprotein

Introduction

Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) is an enveloped, positive-sense, single-stranded RNA virus that causes coronavirus disease 2019 (COVID-19). Virus particles include the RNA genetic material and structural proteins needed for invasion of host cells. Once inside the cell the infecting RNA is used to encode structural proteins that make up virus particles, nonstructural proteins that direct virus assembly, transcription, replication and host control and accessory proteins whose function has not been determined.~ The structural proteins of SARS-CoV-2 include the envelope protein (E), spike or surface glycoprotein (S), membrane protein (M) and the nucleocapsid protein (N). The spike glycoprotein is found on the outside of the virus particle and gives coronavirus viruses their crown-like appearance. This glycoprotein mediates attachment of the virus particle and entry into the host cell. S protein is an important target for vaccine development, antibody therapies and diagnostic antigen-based tests.

Alternative Names

S Mutation (ΔH69/V70, ΔY144, E484K, N501Y, A570D, P681H, T716I, S982A, D1118H, D614G); structural protein; spike protein; Surface glycoprotein; SARS2; SARS-CoV-2; Reporter Virus Particles; RVPs

Gene ID

[43740568](#)

UniProt ID

[P0DTC2](#)