

Product Information

MemDX™ Membrane Protein Mouse C2cd2l (C2 calcium-dependent domain containing 2-like) Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX2042K**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse C2cd2l membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Mouse

Target Protein

C2cd2l

Protein Length

Full Length

Protein Class

Transport

TMD

1

Sequence

MDPDWGQRDVGWAALLVLFAASLITVLGWMLQYARGLWLSRADGGGRDSRPASAAEPGGSLRELGVWRSLLRLRATRTSTPEEAC

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

C2cd2l

Full Name

C2 calcium-dependent domain containing 2-like

Alternative Names

C2cd2l; Tmem24; 1300006O23Rik; phospholipid transfer protein C2CD2L; C2 domain-containing protein 2-like; transmembrane protein 24; C2 calcium-dependent domain containing 2-like

Gene ID

[71764](#)

UniProt ID

[Q80X80](#)