

Product Information

MemDX™ Membrane Protein Mouse Cmtm8 (CKLF-like MARVEL transmembrane domain containing 8) Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX2140K**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse Cmtm8 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Mouse

Target Protein

Cmtm8

Protein Length

Full Length

Protein Class

Cytokine

TMD

4

Sequence

MEEQQRARSHTVTTTTSSFAENFSTTSSSFAYDREFLRTPPGLLIVAEIVLGLLVWTLIAGTEYFRVPAFGWVMFVAVFYWVLTVFFL

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

Cmtm8

Full Name

CKLF-like MARVEL transmembrane domain containing 8

Alternative Names

Cmtm8; Cklfs; Cklfsf8; AA408515; 2700018N07Rik; CKLF-like MARVEL transmembrane domain-containing protein 8; chemokine-like factor super family 8; chemokine-like factor superfamily member 8; CKLF-like MARVEL transmembrane domain containing 8

Gene ID

[70031](#)

UniProt ID

[Q9CZR4](#)