

Product Information

MemDX™ Membrane Protein Mouse Mpv17 (MpV17 mitochondrial inner membrane protein)

Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX2149K**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse Mpv17 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Mouse

Target Protein

Mpv17

Protein Length

Full Length

Protein Class

Receptor

TMD

4

Sequence

MALWRAYQRALAAHPWKVQVLTAGSLMGVGDMISQQQLVERRGLQQHQAGRTLTMVSLGCGFVGPVVGWYKVL D H L I P G T T K V

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

Mpv17

Full Name

MpV17 mitochondrial inner membrane protein

Alternative Names

Mpv17; Tg.Mpv; Tg.Mpv17; protein Mpv17; MpV17 mitochondrial inner membrane protein

Gene ID

[17527](#)

UniProt ID

[P19258](#)