

Product Information

MemDX™ Membrane Protein Mouse Praf2 (PRA1 domain family 2) Expressed *in vitro* *E.coli* expression system, Full Length

Cat. No.: **MPX2152K**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse Praf2 membrane protein expressed *in vitro* *E.coli* expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Mouse

Target Protein

Praf2

Protein Length

Full Length

Protein Class

Transport

TMD

4

Sequence

MSEVRLPPLRALDDFVLGSARLAAPDPGDPQRWCHRVINNLLYYQTNYLLCFGISLALAGYIRPLHTLLSALVVVVALGVLVWAAET

Product Description

Expression Systems

in vitro *E.coli* expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

Praf2

Full Name

PRA1 domain family 2

Alternative Names

Praf2; JM4; Sfc2; Sfc20; C78013; DXImx39; AA986329; DXImx39e; PRA1 family protein 2; PRA1 domain family 2

Gene ID

[54637](#)

UniProt ID

[Q9JIG8](#)