

Product Information

MemDX™ Membrane Protein Mouse Tmem140 (Transmembrane protein 140) Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX2161K**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse Tmem140 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Mouse

Target Protein

Tmem140

Protein Length

Full Length

Protein Class

Receptor

TMD

4

Sequence

MAFSRLWRNNHLPFVGIMILLAAALCLMFYALLWKAGNLADLPSLRIGFYNFCLWKEDMGSLACYNFPELDVLGIPQVGLALARLGV

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

Tmem140

Full Name

Transmembrane protein 140

Alternative Names

Tmem140; 1110007F12Rik; Transmembrane protein 140

Gene ID

[68487](#)

UniProt ID

[Q8BGY5](#)