

Product Information

MemDX™ Membrane Protein Human RABAC1 (Rab acceptor 1) Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX2163K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human RABAC1 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

RABAC1

Protein Length

Full Length

Protein Class

Receptor

TMD

4

Sequence

MAAQKDQQKDAEAEGLSGTTLPLKIPSGAGREWLERRRATIRPWSTFVDQQRFSRPRNLGELCQRLVRNVEYYQSNYVFVFLGL

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

RABAC1

Full Name

Rab acceptor 1

Alternative Names

RABAC1; PRA1; YIP3; PRAF1; prenylated Rab acceptor protein 1; PRA1 domain family 1; PRA1 family protein 1; Rab acceptor 1 (prenylated); prenylated Rab acceptor 1; Rab acceptor 1

Gene ID

[10567](#)

UniProt ID

[Q9UI14](#)