

Product Information

MemDX™ Membrane Protein Human TP53I11 (Tumor protein p53 inducible protein 11)

Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX2169K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human TP53I11 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TP53I11

Protein Length

Full Length

Protein Class

Receptor

TMD

4

Sequence

MAAKQPPPLMKKHSQTDLVSRLKTRKILGVGGEDDDGEVHRSKISQVLGNEIKFTIREPLGLRVWQFVSAVLFSGIAIMALAFPDQLY

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

TP53I11

Full Name

Tumor protein p53 inducible protein 11

Alternative Names

TP53I11; PIG11; tumor protein p53-inducible protein 11; p53-induced gene 11 protein; Tumor protein p53 inducible protein 11

Gene ID

[9537](#)

UniProt ID

[O14683](#)