

Product Information

MemDX™ Membrane Protein Mouse Tex261 (Testis expressed gene 261) Expressed *in vitro* *E.coli* expression system, Full Length

Cat. No.: **MPX2175K**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse Tex261 membrane protein expressed *in vitro* *E.coli* expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Mouse

Target Protein

Tex261

Protein Length

Full Length

Protein Class

Receptor

TMD

5

Sequence

MWFMYYVLSWLSLFIQVAFITLAVAAGLYYLAELIEEYTVATSRIIKYMIWFSTAVLIGLYVFERFPTSMIGVGLFTNLVYFGLLQTFPFIM

Product Description

Expression Systems

in vitro *E.coli* expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

Tex261

Full Name

Testis expressed gene 261

Alternative Names

Tex261; TEG-261; AA409339; AI480706; AL033351; 3110001O07Rik; protein TEX261; testis-expressed protein 261; testis-expressed sequence 261 protein; Testis expressed gene 261

Gene ID

[21766](#)

UniProt ID

[Q62302](#)