

Product Information

MemDX™ Membrane Protein Human NKAIN3 (Sodium/potassium transporting ATPase interacting 3) Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX2176K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human NKAIN3 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

NKAIN3

Protein Length

Full Length

Protein Class

Receptor

TMD

4

Sequence

MGCCTGRCSLICLCALQLVSALERQIFDFLGFWAPILGNFLHIIVVILGLFGTIQYRPRYIMVYTVWTALWVTWNVFIICFYLEVGGLS

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

NKAIN3

Full Name

Sodium/potassium transporting ATPase interacting 3

Introduction

NKAIN3 is a member of a family of mammalian proteins (see NKAIN1; MIM 612871) with similarity to Drosophila Nkain.

Alternative Names

NKAIN3; FAM77D; UG0898H09; NKAIN3-IT1; sodium/potassium-transporting ATPase subunit beta-1-interacting protein 3; NKAIN3 intronic transcript; Na(+)/K(+)-transporting ATPase subunit beta-1-interacting protein 3; Na+/K+ transporting ATPase interacting 3; family with sequence similarity 77, member D; Sodium/potassium transporting ATPase interacting 3

Gene ID

[286183](#)

UniProt ID

[Q8N8D7](#)