

Product Information

MemDX™ Membrane Protein Mouse Tm4sf4 (Transmembrane 4 superfamily member 4)

Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX2186K**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse Tm4sf4 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Mouse

Target Protein

Tm4sf4

Protein Length

Full Length

Protein Class

Receptor

TMD

4

Sequence

MCTGGCARCLGGTLIPLAVFGLLANILLFFPGGKVVNDKSHLSDEVWYFGGILGSGVLMIFPALVFLGLQNNDCCGCCGNEGCGKR

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

Tm4sf4

Full Name

Transmembrane 4 superfamily member 4

Alternative Names

Tm4sf4; il-; Iltmp; transmembrane 4 L6 family member 4; Transmembrane 4 superfamily member 4

Gene ID

[229302](#)

UniProt ID

[Q91XD3](#)