

Product Information

MemDX™ Membrane Protein Human SFT2D3 (SFT2 domain containing 3) Expressed *in vitro* *E.coli* expression system, Full Length

Cat. No.: **MPX2212K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human SFT2D3 membrane protein expressed *in vitro* *E.coli* expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

SFT2D3

Protein Length

Full Length

Protein Class

Transport

TMD

4

Sequence

MADLHRQLQEYLAQGKAGGPAAAEPLLAEEKAEPEGDRPAEEWLGRAGLRWTWARSPAESAAAGLTCLPSVTRGQRLAAGGGC

Product Description

Expression Systems

in vitro *E.coli* expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

SFT2D3

Full Name

SFT2 domain containing 3

Alternative Names

SFT2D3; vesicle transport protein SFT2C; SFT2 domain-containing protein 3; SFT2 domain containing 3

Gene ID

[84826](#)

UniProt ID

[Q587I9](#)