

Product Information

MemDX™ Membrane Protein Mouse Lhfpl3 (Lipoma HMGIC fusion partner-like 3) Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX2227K**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse Lhfpl3 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Mouse

Target Protein

Lhfpl3

Protein Length

Full Length

Protein Class

Receptor

TMD

4

Sequence

MLPAQEAAKLYHTNYVRNSRAIGVLWAIFTICFAINVVCFIQPYWIGDGVDPQAGYFGLFHYCIGNGFSRELTGRGSFTDFSTLPSC

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

Lhfp13

Full Name

Lipoma HMGIC fusion partner-like 3

Alternative Names

Lhfp13; A930031L14Rik; lipoma HMGIC fusion partner-like 3 protein; Lipoma HMGIC fusion partner-like 3

Gene ID

[269629](#)

UniProt ID

[Q9CTN8](#)