

Product Information

MemDX™ Membrane Protein Human CLDN15 (Claudin 15) Expressed *in vitro* *E.coli* expression system, Full Length

Cat. No.: **MPX2245K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human CLDN15 membrane protein expressed *in vitro* *E.coli* expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

CLDN15

Protein Length

Full Length

Protein Class

Transport

TMD

4

Sequence

MSMAVETFGFFMATVGLLMLGVTL PNSYWRVSTVHGNVITNTIFENLWFSCATDSLGVYNCWEFPSMLALSGYIQACRALMITAIL

Product Description

Expression Systems

in vitro *E.coli* expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

CLDN15

Full Name

Claudin 15

Introduction

This gene encodes a member of the claudin family. Claudins are integral membrane proteins and components of tight junction strands. Tight junction strands serve as a physical barrier to prevent solutes and water from passing freely through the paracellular space between epithelial or endothelial cell sheets, and also play critical roles in maintaining cell polarity and signal transductions. Alternatively spliced transcript variants encoding the same protein have been found for this gene.

Alternative Names

CLDN15; claudin-15; Claudin 15

Gene ID

[24146](#)

UniProt ID

[P56746](#)