

Product Information

MemDX™ Membrane Protein Human SCAMP5 (Secretory carrier membrane protein 5)

Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX2269K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human SCAMP5 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

SCAMP5

Protein Length

Full Length

Protein Class

Transport

TMD

4

Sequence

MAEKVNNFPPLPKFIPLKPCFYQDFEADIPPQHVSMTKRLYYLWMLNSVTLAVNLVGCLAWLIGGGGATNFGLAFLWLILFTPCSYV

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

SCAMP5

Full Name

Secretory carrier membrane protein 5

Alternative Names

SCAMP5; secretory carrier-associated membrane protein 5; hSCAMP5; Secretory carrier membrane protein 5

Gene ID

[192683](#)

UniProt ID

[Q8TAC9](#)