

Product Information

MemDX™ Membrane Protein Human STEAP1B (STEAP family member 1B) Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX2288K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human STEAP1B membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

STEAP1B

Protein Length

Full Length

Protein Class

Receptor

TMD

4

Sequence

MESRKDITNQEEIWKMKPRRNLEDNDYLQTAHADEFDCPSELQHAQELFPQWHLPIKIAAVMASLTFLYTLLREVIHPLATSHQQYF

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

STEAP1B

Full Name

STEAP family member 1B

Alternative Names

STEAP1B; STEAP family protein MGC87042; STEAP family member 1B

Gene ID

[256227](#)

UniProt ID

[Q6NZ63](#)