

Product Information

MemDX™ Membrane Protein Human TMEM41A (Transmembrane protein 41A) Expressed *in vitro* E.coli expression system, Full Length of Mature Protein

Cat. No.: **MPX2293K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human TMEM41A membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TMEM41A

Protein Length

Full Length of Mature Protein

Protein Class

Receptor

TMD

5

Sequence

LYLLSTRLPRGRRLGSTEEAGGRSLWFPSDLAELRELSEVLREYRKEHQAYVLLFCGAYLYKQGFAIPGSSFLNVLGALFGPWLC

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TMEM41A

Full Name

Transmembrane protein 41A

Alternative Names

TMEM41A; 2900010K02Rik; Transmembrane protein 41A

Gene ID

[90407](#)

UniProt ID

[Q96HV5](#)