

Product Information

MemDX™ Membrane Protein Mouse Tmem150a (Transmembrane protein 150A) Expressed *in vitro* E.coli expression system, Full Length of Mature Protein

Cat. No.: **MPX2298K**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse Tmem150a membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Mouse

Target Protein

Tmem150a

Protein Length

Full Length of Mature Protein

Protein Class

Receptor

TMD

6

Sequence

MAVMNRHVCPCVENWSYNESCSPDPAEQGGPKSCCTLDDVPLISKCGTYPPECLFSLIGNMGAVMVALICLLRYGQLLEQSRHSW

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

Tmem150a

Full Name

Transmembrane protein 150A

Alternative Names

Tmem150a; Tmem150; BC014685; fasting-inducible integral membrane protein TM6P1; transmembrane protein 150; Transmembrane protein 150A

Gene ID

[232086](#)

UniProt ID

[Q91WN2](#)