

Product Information

MemDX™ Membrane Protein Mouse Cd37 (CD37 antigen) Expressed *in vitro* *E.coli* expression system, Full Length

Cat. No.: **MPX2378K**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse Cd37 membrane protein expressed *in vitro* *E.coli* expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Mouse

Target Protein

Cd37

Protein Length

Full Length

Protein Class

Receptor

TMD

4

Sequence

MSAQESCLSLIKYFLFVNLFVGLGGLIFCFGTWILDKTSFVSFVGLSFVPLQTWSKVLAVSGVLTMALALLGCVGALKELRCLLGL

Product Description

Expression Systems

in vitro *E.coli* expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

Cd37

Full Name

CD37 antigen

Alternative Names

Cd37; Tsp; Tspan26; leukocyte antigen CD37; CD37 antigen

Gene ID

[12493](#)

UniProt ID

[Q61470](#)