

Product Information

MemDX™ Membrane Protein Human ZDHHC18 (Zinc finger DHHC-type palmitoyltransferase

18) Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX2513K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human ZDHHC18 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

ZDHHC18

Protein Length

Full Length

Protein Class

Transferase

TMD

4

Sequence

MKDCEYQQISPGAAPLPASPGARRPGPAASPTPGPGPAPPAAPAPPRWSSSGSGSGSGSLGRRPRRKWEVFPGRNRFYCGO

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

ZDHHC18

Full Name

Zinc finger DHHC-type palmitoyltransferase 18

Alternative Names

ZDHHC18; DHHC18; DHHC-18; palmitoyltransferase ZDHHC18; DHHC domain-containing cysteine-rich protein 18; zinc finger DHHC domain-containing protein 18; zinc finger DHHC-type containing 18; zinc finger, DHHC domain containing 18; Zinc finger DHHC-type palmitoyltransferase 18

Gene ID

[84243](#)

UniProt ID

[Q9NUE0](#)