

Product Information

MemDX™ Membrane Protein Human ZDHHC13 (Zinc finger DHHC-type palmitoyltransferase

13) Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX2832K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human ZDHHC13 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

ZDHHC13

Protein Length

Full Length

Protein Class

Transferase

TMD

6

Sequence

MEGPGLGSQCRNHSHGPHPPGFGRYGICAHENKELANAREALPLIEDSSNCDIVKATQYGIFERCKELVEAGYDVRQPDKENVSL

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

ZDHHC13

Full Name

Zinc finger DHHC-type palmitoyltransferase 13

Alternative Names

ZDHHC13; HIP14L; HIP3RP; palmitoyltransferase ZDHHC13; DHHC-13; HIP14-related protein; huntingtin interacting protein HIP3RP; huntingtin-interacting protein 14-related protein; probable palmitoyltransferase ZDHHC13; putative MAPK-activating protein PM03; putative NF-kappa-B-activating protein 209; zinc finger DHHC domain-containing protein 13; zinc finger DHHC-type containing 13; zinc finger, DHHC domain containing 13; Zinc finger DHHC-type palmitoyltransferase 13

Gene ID

[54503](#)

UniProt ID

[Q8IUH4](#)