

Product Information

MemDX™ Membrane Protein Mouse Zdhhc17 (Zinc finger, DHHC domain containing 17)

Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX2839K**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse Zdhhc17 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Mouse

Target Protein

Zdhhc17

Protein Length

Full Length

Protein Class

Transferase

TMD

6

Sequence

MQREEGFNTKMADGPDEYETETGCVPLLHP EEIKPQSHYNHGYGEPLGRKTHIDDYSTWDIVKATQYGIYERCRELVEAGYDVRQ

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

Zdhhc17

Full Name

Zinc finger, DHHC domain containing 17

Alternative Names

Zdhhc17; Hip; Hip14; DHHC-17; BB187739; A230053P19Rik; D130071N24Rik; palmitoyltransferase ZDHHC17; DHHC domain-containing cysteine-rich protein 17; acyltransferase ZDHHC17; acyltransferase ZDHHC7; huntingtin interacting protein 14; zinc finger DHHC domain-containing protein 17; Zinc finger, DHHC domain containing 17

Gene ID

[320150](#)

UniProt ID

[Q80TN5](#)