

Product Information

MemDX™ Membrane Protein Human TMEM201 (Transmembrane protein 201) Expressed *in vitro* *E.coli* expression system, Full Length

Cat. No.: **MPX2849K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human TMEM201 membrane protein expressed *in vitro* *E.coli* expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TMEM201

Protein Length

Full Length

Protein Class

Receptor

TMD

5

Sequence

MEGVSALLARCPTAGLAGGLGVTACAAAGVLLYRIARRMKPTHTMVNCWFCNQDTLVPYGNRNCWDCPHCEQYNGFQENG DYN

Product Description

Expression Systems

in vitro *E.coli* expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TMEM201

Full Name

Transmembrane protein 201

Alternative Names

TMEM201; Ima1; NET5; SAMP1; RP13-15M17.2; spindle-associated membrane protein 1; Transmembrane protein 201

Gene ID

[199953](#)

UniProt ID

[Q5SNT2](#)