

Product Information

MemDX™ Membrane Protein Human NDC1 (NDC1 transmembrane nucleoporin) Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX2852K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human NDC1 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

NDC1

Protein Length

Full Length

Protein Class

Transport

TMD

6

Sequence

MATAVSRPCAGRSRDILWRVLGWRIVASIVWSVLFLPICTTVFIIFSRIDLFHPIQWLSDSFSDLYSSYVIFYFLLLSVVIISIFNVEFYA

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

NDC1

Full Name

NDC1 transmembrane nucleoporin

Alternative Names

NDC1; NET3; TMEM48; nucleoporin NDC1; nuclear division cycle 1 homolog; transmembrane protein 48; NDC1 transmembrane nucleoporin

Gene ID

[55706](#)

UniProt ID

[Q9BTX1](#)