

Product Information

MemDX™ Membrane Protein Human VN1R4 (Vomeronasal 1 receptor 4) Expressed *in vitro* *E.coli* expression system, Full Length

Cat. No.: **MPX2931K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human VN1R4 membrane protein expressed *in vitro* *E.coli* expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

VN1R4

Protein Length

Full Length

Protein Class

GPCR

TMD

7

Sequence

MASRYVAVGMILSQTVVGVLGFSVLLHYLSFYCTGCRLRSTDLIVKHLIVANFLALRCKGVPQTMAAFGVRYFLNALGCKLVFYLFH

Product Description

Expression Systems

in vitro *E.coli* expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

VN1R4

Full Name

Vomeronasal 1 receptor 4

Alternative Names

VN1R4; V1RL4; CTD-2245F17.9; G-protein coupled receptor GPCR27; V1R-like 4; V1r-like receptor 4; hGPCR27; pheromone receptor; Vomeronasal 1 receptor 4

Gene ID

[317703](#)

UniProt ID

[Q7Z5H5](#)