

Product Information

MemDX™ Membrane Protein Human OR4A16 (Olfactory receptor family 4 subfamily A member 16) Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX3333K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human OR4A16 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

OR4A16

Protein Length

Full Length

Protein Class

GPCR

TMD

7

Sequence

MRPSSNVTEFVLLGLTQDPDVKKTLFVMFLLIYIVTMVGNLLIWVTTIGSPSLGSLMYFFLAYLSLMDAIYSTAMSPKLMIDLLCDKIAIS

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

OR4A16

Full Name

Olfactory receptor family 4 subfamily A member 16

Introduction

Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms.

Alternative Names

OR4A16; OR4A16Q; OR11-117; olfactory receptor 4A16; olfactory receptor OR11-117; Olfactory receptor family 4 subfamily A member 16

Gene ID

[81327](#)

UniProt ID

[Q8NH70](#)