

Product Information

MemDX™ Membrane Protein Human MCHR2 (Melanin concentrating hormone receptor 2)

Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX3381K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human MCHR2 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

MCHR2

Protein Length

Full Length

Protein Class

GPCR

TMD

7

Sequence

MNPFHASCWNTSAELLNKS WNKEFAYQTASVVDTVILPSMIGIICSTGLVGNILIVFTIIRSRKKTVPDIYICNLAVADLVHIVGMPFLIH

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

MCHR2

Full Name

Melanin concentrating hormone receptor 2

Alternative Names

MCHR2; SLT; MCH2; MCH2R; GPR145; GPRv17; MCH-2R; MCH-R2; MCHR-2; G protein-coupled receptor slt; G-protein coupled receptor 145; MCH receptor 2; melanin-concentrating hormone 2; Melanin concentrating hormone receptor 2

Gene ID

[84539](#)

UniProt ID

[Q969V1](#)