

Product Information

MemDX™ Membrane Protein Human MRGPRF (MAS related GPR family member F)

Expressed *in vitro* *E.coli* expression system, Full Length

Cat. No.: **MPX3390K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human MRGPRF membrane protein expressed *in vitro* *E.coli* expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

MRGPRF

Protein Length

Full Length

Protein Class

GPCR

TMD

7

Sequence

MAGNCSWEAHPGNRNKMCPGLSEAPELYSRGFLTIEQIAMLPPPAVMNYIFLLLCLCGLVGNGLVLWFFGFSIKRNPFSIYFLHLAS

Product Description

Expression Systems

in vitro *E.coli* expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

MRGPRF

Full Name

MAS related GPR family member F

Alternative Names

MRGPRF; RTA; MRGF; GPR140; GPR168; mas-related G-protein coupled receptor member F; G protein-coupled receptor MrgF; G-protein coupled receptor 140; G-protein coupled receptor 168; MAS-related GPR, member F; Mas-related G protein-coupled receptor F; mas-related G protein-coupled MRGF; mas-related gene F protein; MAS related GPR family member F

Gene ID

[116535](#)

UniProt ID

[Q96AM1](#)