

Product Information

MemDX™ Membrane Protein Human PAQR6 (Progestin and adipoQ receptor family member 6) Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX3394K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human PAQR6 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

PAQR6

Protein Length

Full Length

Protein Class

Receptor

TMD

7

Sequence

MLSLKLPQLLQVHQVPRVFWEDGIMSGYRRPTSSALDCVLSSFQMTNETVNIWTHFLPTWYFLWRLLALAGGPGFRAEPYHWPLL

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

PAQR6

Full Name

Progesterin and adipoQ receptor family member 6

Alternative Names

PAQR6; PRdelta; membrane progesterin receptor delta; mPR delta; membrane progesterone P4 receptor delta; membrane progesterone receptor delta; progesterone and adipoQ receptor family member 6; progesterin and adipoQ receptor family member VI variant 3; progesterin and adipoQ receptor family member VI variant 4; progesterin and adipoQ receptor family member VI variant 5; Progesterin and adipoQ receptor family member 6

Gene ID

[79957](#)

UniProt ID

[Q6TCH4](#)