

Product Information

MemDX™ Membrane Protein Human FPR3 (Formyl peptide receptor 3) Expressed *in vitro* *E.coli* expression system, Full Length

Cat. No.: **MPX3429K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human FPR3 membrane protein expressed *in vitro* *E.coli* expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

FPR3

Protein Length

Full Length

Protein Class

GPCR

TMD

7

Sequence

METNFSIPLNETEEVLPEPAGHTVLWIFSLLVHGVTFVFGVLGNGLVIWVAGFRMTRTVNTICYLNALADFSFSAILPFRMVSVAMR

Product Description

Expression Systems

in vitro *E.coli* expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

FPR3

Full Name

Formyl peptide receptor 3

Alternative Names

FPR3; FMLPY; FPRH1; FPRH2; FPRL2; RMLP-R-I; FMLP-R-II; FML2_HUMAN; N-formyl peptide receptor 3; FMLP-related receptor II; formyl peptide receptor-like 2; Formyl peptide receptor 3

Gene ID

[2359](#)

UniProt ID

[P25089](#)