

Product Information

MemDX™ Membrane Protein Human GPR37L1 (G protein-coupled receptor 37 like 1)

Expressed *in vitro* *E.coli* expression system, Full Length of Mature Protein

Cat. No.: **MPX3588K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human GPR37L1 membrane protein expressed *in vitro* *E.coli* expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

GPR37L1

Protein Length

Full Length of Mature Protein

Protein Class

GPCR

TMD

7

Sequence

APLHLGRHRAETQEQQSRSKRGTEDEEAKGVQQYVPEEWAHEYPRPIHPAGLQPTKPLVATSPNPGKDGGTPDSGQELRGNLTGA

Product Description

Expression Systems

in vitro *E.coli* expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

GPR37L1

Full Name

G protein-coupled receptor 37 like 1

Alternative Names

GPR37L1; ETBRLP2; ETBR-LP-2; ET(B)R-LP-2; endothelin B receptor-like protein 2; endothelin type b receptor-like protein 2; prosaposin receptor GPR37L1; G protein-coupled receptor 37 like 1

Gene ID

[9283](#)

UniProt ID

[O60883](#)