

Product Information

MemDX™ Membrane Protein Human CHRM5 (Cholinergic receptor muscarinic 5) Expressed *in vitro* E.coli expression system, Full Length

Cat. No.: **MPX3623K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human CHRM5 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

CHRM5

Protein Length

Full Length

Protein Class

GPCR

TMD

7

Sequence

MEGDSYHNATTVNGTPVNHQPLERHRLWEVITIAAVTAVVSLITIVGNVLVMISFKVNSQLKTVNNYYLLSLACADLIIGIFSMNLYTTY

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

CHRM5

Full Name

Cholinergic receptor muscarinic 5

Introduction

The muscarinic cholinergic receptors belong to a larger family of G protein-coupled receptors. The functional diversity of these receptors is defined by the binding of acetylcholine and includes cellular responses such as adenylate cyclase inhibition, phosphoinositide degeneration, and potassium channel mediation. Muscarinic receptors influence many effects of acetylcholine in the central and peripheral nervous system. The clinical implications of this receptor are unknown; however, stimulation of this receptor is known to increase cyclic AMP levels.

Alternative Names

CHRM5; HM5; muscarinic acetylcholine receptor M5; acetylcholine receptor, muscarinic 5; Cholinergic receptor muscarinic 5

Gene ID

[1133](#)

UniProt ID

[P08912](#)