

Product Information

MemDX™ Membrane Protein Human GPR107 (G protein-coupled receptor 107) Expressed *in vitro* E.coli expression system, Full Length of Mature Protein

Cat. No.: **MPX3634K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human GPR107 membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

GPR107

Protein Length

Full Length of Mature Protein

Protein Class

GPCR

TMD

7

Sequence

RVHHLALKDDVRHKVHLNTFGFFKGYMNVVSSLSLNEPEDKDVITIGFSLDRTKNDGFSSYLDEDVNYCILKKQSVSVTLLILDISR

Product Description

Expression Systems

in vitro E.coli expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

GPR107

Full Name

G protein-coupled receptor 107

Alternative Names

GPR107; GCDRP; LUSTR1; bA138E2.2; protein GPR107; lung seven transmembrane receptor 1; G protein-coupled receptor 107

Gene ID

[57720](#)

UniProt ID

[Q5VW38](#)